
● CHROMagar™ Salmonella

For detection and isolation of *Salmonella*

For detection and isolation of *Salmonella* species, including *S.typhi* and *S.paratyphi* in clinical specimens

Background

Infections caused by *Salmonella* spp, including *Salmonella typhi*, remain a major worldwide health problem:

- In the US, *Salmonella* has an incidence rate of 16.2 cases per 100,000 (CDC estimation, 2008).
 - In Europe, it is reported as the first cause of collective toxi-infections. (2007 EFSA report)
 - In developing countries, *Salmonella typhi* and *paratyphi* are commonly encountered with an estimated annual incidence of about 17 million cases. (2007 EFSA report)
- Moreover, according to a recent WHO report, *Salmonella* infections are responsible for 2 million deaths per year from diarrhoea.

Mainly due to contamination in the food chain and/or during food-production processes, *Salmonella* commonly induces enteric illness whose major symptoms are abdominal cramps, diarrhea, nausea, vomiting. More severe cases, for instance typhoid cases or infections in immuno-depressed patients, can lead to body dehydration with renal failure or bacteraemia.

Medium Performance

1 EASY READING

intense mauve colony colours for better identification and partial inhibition of *E.coli* and coliforms.

2 GREATER SPECIFICITY / LESS WORKLOAD

conventional media for the detection of *Salmonella* by H₂S character have very poor specificity resulting in numerous false positives (*Citrobacter*, *Proteus*, etc.) among the rare, real positive *Salmonella*. The workload for unnecessary examination of suspect colonies is so heavy that real positive *Salmonella* colonies might often be overlooked in routine testing. Because of their poor specificity, conventional media require a tedious examination of at least 10 colonies per suspected sample. On the contrary, CHROMagar™ Salmonella eliminates most of those false positives and allows technicians to focus on the real contaminated samples.

3 HIGH SENSITIVITY AND SPECIFICITY

leading to a higher detection rate of *Salmonella*

Sensitivity: 100%*

Specificity: 89%* compared to 78% with Hektoen Agar.

* Specificity and sensitivity from scientific study: "Comparison of CHROMagar Salmonella medium and Hektoen Enteric Agar for isolation of *Salmonella* from stool samples." Gaillot O. et al. 1999. Journal of Clinical Microbiology, 37 : 762-765

4 DRAMATIC REDUCTION OF THE WORKLOAD

Number of useless confirmatory tests is minimized since there is no need of duplicating them.

Medium Description

Powder Base	Total	34.9 g/L
	Agar	15.0
	Peptone & Yeast extract	7.0
	Chromogenic and selective mix	12.9
	Storage at 15/30°C - pH: 7.6 ± 0.2	
	Shelf Life	2 years

Usual Samples	- Syndrom typhoid → stool or blood samples - Gastro enteritis → stool samples
Procedure	Direct Streaking. Incubation 24h, 37°C. Aerobic condition.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com



Plate Reading

- *Salmonella* including *S.typhi* → mauve
- Other bacteria → blue, colourless or inhibited.

Quality Control Strains

<i>S. enteritidis</i> ATCC® 13076	mauve
<i>S. typhimurium</i> ATCC® 13311	mauve
<i>E. coli</i> ATCC® 25922	metallic blue, small
<i>C. freundii</i> ATCC® 8090	metallic blue
<i>C. albicans</i> ATCC® 60193	inhibited
<i>S. aureus</i> ATCC® 25923	inhibited

ATCC® is a registered trademark of the American Type Culture Collection


Order References

Please use these product references when contacting your local distributor:

- 1000 ml pack SA130
- 5000 ml pack SA132
- 25 L pack SA133-25
- Bulk on request

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com
Find your nearest distributor on
www.CHROMagar.com/contact

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